

ANTIBODY SPECIFICATION

CAT. No. 11121

Products : Recombinant Proteins for HCG Capture

Descriptions: Recombinant Proteins for HCG Capture

Clone Name: 3B9

Isotype: IgG2a

Specificity: Hepatitis B Virus(HBV) specific.

Immunogen: Native Hepatitis B Virus.

Host animal: Mouse. Hybridoma from Sp2/0 myeloma cell and lymphocytes in Balb/c mice.

Format: Purified immunoglobulin, Liquid.

Purification: 95% or above purified, Protein G sepharose chromatography.

Buffer: PBS (pH 7.4) with 0.1% Sodium azide.

Applications: Enzyme Linked Immunosorbent Assay (ELISA).
Lateral Flow Immunoassay (LFIA).
Optimum working dilutions must be determined by user.

Storage: Store at 2~8°C,
Antibody solutions can be stored at -20°C but avoid repeated freezing and thawing.

In Vitro use only.
Research use only