

ANTIBODY SPECIFICATION

CAT. No. 11122

Products : Anti-HBsAg

Descriptions: Monoclonal antibody for the detection of HBsAg

Clone Name: 6F52

Isotype: IgG1

Specificity: Hepatitis B Virus(HBV) specific.

Immunogen: Native Hepatitis B Virus.

Host animal: Mouse. Hybridoma from Sp2/0 myeloma cell and lymphocytes in Balb/c mice.

Format: Purified immunoglobulin, Liquid.

Purification: 95% or above purified, Protein G sepharose chromatography.

Buffer: PBS (pH 7.0) with 0.1% Sodium azide.

Applications: Enzyme Linked Immunosorbent Assay (ELISA).
Lateral Flow Immunoassay (LFIA).
Optimum working dilutions must be determined by user.

Storage: Store at 2~8°C,
Antibody solutions can be stored at -20°C but avoid repeated freezing and thawing.

In Vitro use only.
Research use only

878 Sunhwan-daero, Chuncheon-si, Gangwon-do 24399, Republic of Korea

Tel: +82 (0)33 244 0100 Fax: +82 (0)33 244 4634

E-mail: median@mediandx.com