

ANTIBODY SPECIFICATION

CAT. No. 17401

Products : Anti-Hemoglobin A0

Descriptions: Monoclonal Antibody to Hemoglobin.

Clone Name: M230492

Isotype: IgG1

Specificity: Human Hemoglobin specific.

Immunogen: Human Hemoglobin.

Host animal: Mouse. Hybridoma from Sp2/0 myeloma cell and lymphocytes in Balb/c mice.

Format: Purified immunoglobulin, Liquid.

Purification: 95% or above purified, Protein G sepharose chromatography.

Buffer: PBS (pH 7.0) with 0.1% Sodium azide.

Applications: Enzyme Linked Immunosorbent Assay (ELISA).
Lateral Flow Immunoassay (LFIA).
Optimum working dilutions must be determined by user.

Storage: Store at 2 ~ 8°C,
Antibody solutions can be stored at -20°C but avoid repeated freezing and thawing.

In Vitro use only.
Research use only