



## PRODUCT SPECIFICATION

CATALOG NUMBER : 22301

PRODUCT	<b>rec. Malaria Pf MSP 1-1</b> [ <i>Plasmodium falciparum</i> merozoite surface protein]
PRODUCT DESCRIPTION	The recombinant Pf-MSP 1-1 was expressed in <i>E. coli</i> and purified by affinity chromatography to homogeneity. The purified recombinant protein was concentrated and dialyzed with phosphate-buffered saline.
CATALOG No.	22301
EXPRESSED IN	<i>E. coli</i>
PRESERVATIVES	0.05% Sodium Azide
BUFFER	10mM PBS, pH 7.2
BIOHAZARD	Use normal laboratory precautions and procedures.
APPLICATIONS	Rapid Diagnostic Test(RDT), EIA
PURITY	>95% by SDS-PAGE gel
STORAGE	Store at -20°C.

For In Vitro Use Only