

ANTIBODY SPECIFICATION

CAT. No. 9042

Products : Anti-PRRSV Monoclonal Antibody

Descriptions: Monoclonal Antibody to Porcine reproductive and respiratory syndrom virus (PRRSV)

NC protein (European type specific).

Clone Name: 6F34

Isotype: IgG2a

Specificity: Specific to PRRSV NC protein of European type.
No or a little reactive to PRRSV American type (VR-2332 strain).

Immunogen: Purified recombinant NC protein from Lelystad strain (European type).

Host animal: Mouse. Hybridoma from Sp2/0 myeloma cell and lymphocytes in Balb/c mice.

Format: Purified immunoglobulin, Liquid.

Purification: 95% or above purified, Protein G sepharose chromatography.

Buffer: PBS (pH 7.4) with 0.1% Sodium azide.

Application: Fluorescent Antibody (FA) Test,
Enzyme Linked Immunosorbent Assay (ELISA).
Optimum working dilutions must be determined by user.

Storage: Store at 2~8°C

Antibody solutions can be stored at -20°C but avoid repeated freezing and thawing.

In Vitro use only.

Research use only

878 Sunhwan-daero, Chuncheon-si, Gangwon-do 24399, Republic of Korea

Tel: +82 (0)33 244 0100 Fax: +82 (0)33 244 4634

E-mail: median@mediandx.com