

## ANTIBODY SPECIFICATION

**CAT. No. 9162**

**Products : Anti-CCV Monoclonal Antibody**

**Descriptions:** Monoclonal Antibody to Canine Corona Virus (CCV).

**Clone Name:** M32345

**Isotype:** IgG1

**Specificity:** Canine Corona Virus (CCV) specific.  
No cross react to CPV, CDV.

**Immunogen:** Purified CCV virus (ATCC VR-809).

**Host animal:** Mouse. Hybridoma from Sp2/0 myeloma cell and lymphocytes in Balb/c mice.

**Format:** Purified immunoglobulin, Liquid.

**Purification:** 95% or above purified, Protein G sepharose chromatography.

**Buffer:** PBS (pH 7.4) with 0.1% Sodium azide.

**Applications:** Fluorescent Antibody (FA) Test, Immunohistochemical (IHC) staining.  
Lateral Flow Immunoassay (LFIA).  
Optimum working dilutions must be determined by user.

**Storage:** Store at 2~8°C,  
Antibody solutions can be stored at -20°C but avoid repeated freezing and thawing.

*In Vitro use only.  
Research use only*