

ANTIBODY SPECIFICATION

CAT. No. 9165

Products : Anti-CDV Monoclonal Antibody

Descriptions: Monoclonal Antibody to Canine Distemper Virus (CDV).

Clone Name: 9A5

Isotype: IgG2b

Specificity: Canine Distemper Virus (CDV) specific.
No cross react to CPV, CCV.

Immunogen: Purified recombinant H protein of CDV.

Host animal: Mouse. Hybridoma from Sp2/0 myeloma cell and lymphocytes in Balb/c mice.

Format: Purified immunoglobulin, Liquid.

Purification: 95% or above purified, Protein G sepharose chromatography.

Buffer: PBS (pH 7.4) with 0.1% Sodium azide.

Applications: Fluorescent Antibody (FA) Test.
Enzyme Linked Immnosorbent Assay (ELISA).
Optimum working dilutions must be determined by user.

Storage: Store at 2~8°C,
Antibody solutions can be stored at -20°C but avoid repeated freezing and thawing.

In Vitro use only.
Research use only