

## **POLYCLONAL ANTIBODY SPECIFICATION**

**CAT. No. 9901**

**Products : Anti-CSFV Rabbit Poly-Antibody**

**Descriptions:** Purified rabbit antibody immunized with CSFV recombinant E2 protein.

**Specificity:** Reactive to Classical Swine Fever Virus E2 protein.  
Low level cross reactive to BVDV virus in Fluorescent antibody(FA) test.

**Immunogen:** Purified CSFV recombinant E2 protein.

**Host animal:** Serum from Rabbit, Newzealand white.

**Purification:** Protein A affinity chromatography.

**Buffer:** PBS (pH 7.4) with 0.1% Sodium azide.

**Applications:** Fluorescent Antibody (FA) Test,  
Enzyme Linked Immnosorbent Assay (ELISA).  
Optimum working dilutions must be determined by user.

**Storage:** Store at 2~8°C,  
Antibody solutions can be stored at -20°C but avoid repeated freezing and thawing.

*In Vitro use only.*  
*Research use only*

878 Sunhwan-daero, Chuncheon-si, Gangwon-do 24399, Republic of Korea

Tel: +82 (0)33 244 0100 Fax: +82 (0)33 244 4634

E-mail: median@mediandx.com