

## **POLYCLONAL ANTIBODY SPECIFICATION**

**CAT. No. 9903**

**Products : Anti-PRRSV Rabbit Poly-Antibody**

**Descriptions:** Purified rabbit antibody immunized with PRRSV recombinant NC protein.

**Specificity:** NC protein of Pocrine reproductive and respiratory virus (North american strain). Cross reactive to PRRSV European strain.

**Immunogen:** Purified recombinant NC protein of PRRSV strain VR2332.

**Host animal:** Serum from Rabbit, Newzealand white.

**Purification:** Protein A affinity chromatography.

**Buffer:** PBS (pH 7.4) with 0.1% Sodium azide.

**Applications:** Fluorescent Antibody (FA) Test,  
Enzyme Linked Immnosorbent Assay (ELISA).  
Optimum working dilutions must be determined by user.

**Storage:** Store at 2~8°C,  
Antibody solutions can be stored at -20°C but avoid repeated freezing and thawing.

*In Vitro use only.*  
*Research use only*