

VDRG® E.canis Ab Rapid Kit



Ehrlichia canis

Ehrlichia is a rickettsial genus of bacteria belonging to the family Ehrlichiaeae. There are several species of Ehrlichia, but the one that most commonly affects dogs and causes the most severe clinical signs is Ehrlichia canis. This test kit, the diagnostic reagent can detect Ehrlichia Canis antibodies quickly and simply within 10 minutes after dropping the samples.

Introduction


- + **Intended use** : Detection of Ehrlichia canis antibody
- + **Principle** : Immunochromatographic assay
- + **Specimen** : Canine whole blood, plasma or serum
- + **Component**
 - E.canis Ab Rapid device
 - Dilution buffer bottle
 - Capillary tube
 - Anticoagulant tube



Features

- + Clinical Sensitivity : 98.3%(121/123)
- + Clinical Specificity : 100%(305/305)

Test Procedure



Cat No.	Product Name	Quantity
PC-ECA-11	VDRG® E.canis Ab Rapid Kit	10 Tests/Box